0m+

LYRA

digital baby scale

<b>EN</b>	INSTRUCTIONS FOR USE	<b>AL</b>	UDHËZIMET E PËRDORIMIT
<b>ES</b>	INSTRUCCIONES DE USO	<b>BS</b>	UPUTSTVO ZA UPOTREBU
<b>BG</b>	ИНСТРУКЦИИ ЗА УПОТРЕБА	<b>KA</b>	გამოყენების ინსტრუქცია
<b>IT</b>	ISTRUZIONI D'USO	<b>UA</b>	ІНСТРУКЦІЯ З ВИКОРИСТАННЯ
<b>GR</b>	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	<b>HR</b>	UPUTE ZA UPORABU
<b>DE</b>	GEBRAUCHSANWEISUNG	<b>CZ</b>	NÁVOD K POU ITÍ
<b>FR</b>	MODE D'EMPLOI	<b>SK</b>	NÁVOD NA POU ITIE
<b>NL</b>	GEBRUIKSAANWIJZING	<b>SL</b>	NAVODILA ZA UPORABO
<b>FI</b>	KÄYTTÖOHJEET	<b>HU</b>	HASZNÁLATI UTASÍTÁS
<b>RO</b>	INSTRUCȚIUNI DE UTILIZARE	<b>HY</b>	ՕԳՏԱԳՈՐԾՄԱՆ ՑՈՒՑՈՒՄՆԵՐ
<b>PT</b>	INSTRUÇÕES DE USO	<b>ET</b>	KASUTUSJUHEND
<b>MT</b>	ISTRUZZJONIJIET GHALL-UŻU	<b>LT</b>	NAUDOJIMO INSTRUKCIJA
<b>PL</b>	INSTRUKCJA UŻYTKOWANIA	<b>LV</b>	LIETOŠANAS INSTRUKCIJA
<b>SR</b>	УПУТСТВО ЗА УПОТРЕБУ	<b>FA</b>	دستور العمل استفاده
<b>MK</b>	УПАТСТВА ЗА УПОТРЕБА	<b>AR</b>	تعليمات الاستخدام

## INSTRUCTIONS FOR USE

Dear customer,

Thank you for choosing a product of our range. Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

## WARNINGS

- Never leave the child unattended
- Protect the baby scale from knocks, moisture, dust, chemicals and drastic changes in temperature.
- Keep packaging material away from children (risk of suffocation)
- Do not place any objects on the scale when not in use
- Remove any transport locks
- The precision of the scale may be impaired by strong electromagnetic fields (e.g. induction, hobs, mobile phones)

## WARNINGS

### BATTERY

- If your skin or eyes come into contact with battery fluid, rinse the affected area with water and seek medical assistance.
- Small children may swallow and choke on batteries
- If the battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth
- Protect batteries from excessive heat
- Do not throw batteries into a fire
- Use identical or equivalent battery types only
- Do not disassemble, open or crush the batteries

## SPECIFICATIONS

• Equipped with a High Precision "Strain-Gauge" Sensor • Large backlit LCD display • Tare weight function • Max capacity: 20kg • Graduation: 5g • Adjustable kg/lb/oz • Automatic switch off • 2xAAA batteries included

The scale has four anti-slip pads on bottom for safety use.

The whole set includes also a measuring tape for baby growth and 2xAAA batteries.

## IMPORTANT NOTES

### 1. Caution

- Protect your personal scales from impact with hard objects, moisture, dust, chemicals, toiletries, liquid cosmetics, great temperature fluctuation and closeness to sources of heat (open fires, radiators).
- Cleaning: You can clean the scale with a damp cloth and a little washing up liquid, if required. Never immerse the scale in water or rinse it under running water.
- Storage: Do not place and objects on the scale when not in use.
- Strong electromagnetic fields (e.g. Cell phones) may impair the accuracy of the scale.
- The scales can accept a maximum load of 44 lb (20 kg/704 oz). The weight measurement is displayed in 0.20 oz (5 g/0.02 lb) divisions.
- The scale is not intended to be used for commercial purposes.
- Remove transport safety catch.
- Place the scale on a level and solid surface; this is important for ensuring that the scale weighs correctly.

### 2. Operation

#### 2-1 Battery

The scale required 2x AAA batteries. Before using the baby scales, these must be inserted in the battery compartment beneath the scales, ensuring the right polarity.

When "Lo" appears in the display, the battery need to be replaced.

#### 2-2 Weighing

The weight is displayed in steps of 5g from a minimum weight of 50g.

Press the button "TARE" to switch on the scale.

As soon as "0.000 kg" appears on the display, the scale is ready to take measurements.

Now place the baby on the weighing area. Do this very carefully; do not let the baby fall and take particular care with the baby's head.

This baby scale has a "HOLD" function (AUTO HOLD and MECHANICAL HOLD). The AUTO HOLD function establishes the baby's weight on the display as soon as he/she is still enough to achieve an exact measurement value. If the baby is moving around vigorously, his/her weight cannot be recorded.

## CARE AND ATTENTION

- In normal weighing, the scale only receives static pressure. If the scale is used abnormally by impact, or dropped pressure, the glass may break easily and may cause danger.
- Please stand evenly on the scale and do not move.
- Do not use any chemical cleaning agent. Clean your scale with damp cloth, but do not allow water to get inside the scale.
- Put the scale in a cool and dry place. Do not store it in wet environment. As it may destroy the battery and shorten the battery life.
- Put the scale on flat surface; don't put it in an upright position.
- If there is not display after the scale starts, please check whether the battery and battery contacts are well connected. If the scale has not been used for a long period, please check that you need to replace the battery or not.
- The scale is for family weighing use only not for commercial or other purposes.
- Please do not put heavy objects on the scale.

Therefore, try to calm the wriggling baby to achieve a rapid measurement, but do not touch the baby, as this will affect the weight measurement.

As soon as the baby's weight is established on the display, the word "HOLD" stops flashing and the weight is displayed for a further 20 seconds. The scale then switches off automatically.

If the baby is wriggling about too vigorously, you can select the MECHANICAL HOLD function by pressing the "HOLD" button. Within 3 seconds, the word "HOLD" stops flashing and the weight is displayed for a further 20 seconds.

If you do not want to place the baby directly on the weighing surface, but prefer to use a towel, for example, you have to place the towel on the scale BEFORE switching them on. If the towel (and any other items which you want to deduct again during weighing) weighs over 50 g, you can also use the TARE function. Once you have pressed the TARE button, 'TARE' is shown in the display. To reset the scale to 0.000g, press the TARE button after "HOLD" has stopped flashing.

On delivery the scale is set at "kg" units. The button "kg/lb/oz" allows you to switch at any time between kilograms units and pounds(lb) and ounces(oz).

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Kika Group Ltd.